

THE FLOWRAMAN CELL

with BallProbe technology

Wetted Materials

Probe Body	Hastelloy C22
Optical Interface	6mm diameter UV-grade sapphire BallProbe
Sealing Materials	Perfluoroelastomer (Kalrez®)
	Gold option for higher pressure applications

Specifications

Fluid Connectors	Compatible with 1/8 in. Swagelok or Parker A-lok fittings
Body Dimensions	W 35mm H 20mm D 13mm
Fiber Probe Interface	80mm long optical connection compatible with MarqMetrix Fiber BallProbe.
	Also available with 80mm long 0.5 in. (12.7mm) straight tube connection.
Continuous Operating Temperature Range	-20°C to 250°C with perfluoroelastomer seal
	-20°C to 350°C with gold seal
Individual Pressure Test Rating	500 psi (34 bar) with perfluoroelastomer seal
	2500psi (170 bar) with gold seal
Compatible Laser Wavelengths	500-1100nm
Flow Path Volume	200uL

The MarqMetrix® FlowRaman® Cell

is a Raman BallProbe® intended for demanding process flow-through applications. The unibody design places the patented BallProbe spherical optic directly in the flow path, minimizing dead volume and sample handling requirements while taking full advantage of the optical efficiencies.

The MarqMetrix FlowRaman Cell (FRC) is designed for easy integration into continuous flow applications. Plumbing connection are made with standard 1/8 in. tube fitting connectors (e.g. Swagelok, Parker A-lok, Hoke, etc.) Each FRC is individually pressure tested so that you can quickly and confidently integrate powerful Raman measurements into your continuous flow process.



A medium-pressure FRC is available compatible with pressures up to 15,000psi, and temperatures up to 400°C. Please call MarqMetrix to discuss your application requirements.

For more information call:

206-971-3625

info@MarqMetrix.com
MarqMetrix.com

The power of Raman, the simplicity of MarqMetrix

THE FLOWRAMAN CELL

with BallProbe immersion technology



The use of high-performance Hastelloy and sapphire makes the FRC resistant to harsh chemical environments. The FRC has been installed in commercial applications ranging from cryogenic environments to boiling concentrated acids.

Additionally, the curved shape of the sapphire ball facilitates material exchange near the surface, preventing the buildup of materials that interfere with spectral acquisition. The form factor and 'self-cleaning' properties of the design make the FRC an ideal solution for Raman measurements in process flow applications.

The FRC is optimized to integrate with the MarqMetrix Fiber BallProbe, but can be used with other fiber optic probes that have a collimated excitation beam co-linear with the collection path.

Optical Properties

Optimized for use with the MarqMetrix Fiber Probe

Accepts a collimated laser beam up to 5.6mm (0.22 in.) diameter

Made with high purity UV-grade sapphire ball lens aligned along the C-axis, eliminating response variability due to birefringence

Operating Conditions

Suitable for continuous exposure to dilute and concentrated acids (hot & cold), bases and most organic solvents including ethanol, THF, ethyl acetate, acetone, DCM, toluene, pentane and acetonitrile

Avoid prolonged exposure to aqua regia (gold option)


Related Products

Process BallProbe - 0.5 in. diameter, for more demanding environment

Fiber BallProbe - filtered fiber optic interface specifically designed for the MarqMetrix BallProbe and MarqMetrix Raman All-In-One.



High-performance Hastelloy and sapphire makes the FRC impervious to harsh chemical environments.



For more information call:


206-971-3625

info@MarqMetrix.com
MarqMetrix.com

8021. v1.0

MARQMETRIX

U.S. Patent No. 6,977,729 B2

MarqMetrix BallProbe TouchRaman  are registered trademarks of MarqMetrix, Inc.